

QUALITY FOR LIFE



A200 Power Wheelchair

Get in and go!

The A200 is the perfect consumer power wheelchair for anyone who simply wants to get in and go! Its rear-wheel drive makes it easy for first time users to control, and at only 22 ½ narrow inches, its stable, compact design maneuvers nimbly in any indoor environment.

Tool-less seating adjustments and a swing-away joystick allow the A200 to be quickly adapted to the user, while mechanical elevating leg rests easily lift legs for a change in position. Want to visit friends and family? Modular construction and integral tie-down points make the A200 a great travel partner.

The A200. Get in and go.



A200 Power Wheelchair



Contoured seat (option)

A200 Standard Features and Benefits

- Modular construction combined with direct drive for easy setup and service
- Compact design, easily navigates through indoor environment
- Tool-less adjustments, easily adapted to driver
- Integrated tie-down points for increased transportation options¹

HCPCS Codes:²

K0822 K0823



Seat angle easily adjusted

Penny & Giles (P&G) VR2 integral controller

Technical Data	
Maximum speed	4 mph
Battery capacity	U1, 31 Ah
Distance range	approx. 12 miles ³
Maximum charging time	12 hr
Max. load	220 lbs
Weight when empty	145 lbs
Max. obstacle	2"
Total width	23″
Total length without Legrests	31½″
Total length with Legrests	39½″
Seat height (seat set at 0°)	17" or 19"
Seat width	15" – 21½"
Seat depth	15" – 18"
Seat inclination	0°/3°/6° (adjustable without tools)
Armrest height	9½″ – 14″
Back angle (adjustable)	-9°/1°/11°/21° or 0°/10°/20°/30°
Lower leg length	10" – 17½" (add 3½" for contoured seat)

¹ The A200 has been tested to and passed RESNA's WC19 transportation standard for "wheelchairs used as seats in motor vehicles". Otto Bock recommends when ever possible all passengers using public transportation transfer to the bus seat provided.

² The responsibility for accurate coding lies with the patient care facility that selects a product



² The responsibility for accurate coding lies with the patient care facility that selects a product and fits a patient. Otto Bock's coding recommendations are based on our best judgment. These recommendations are open to revision based on additional information or changes in the alpha-numeric system.

³ Your driving range may differ depending on terrain and driving style.

Specifications and prices are subject to change without notice.