

Duomed

The economic solution for venous problems.

medi tip¹:

Due to the material thickness, the fit and the product characteristics, we recommend Duomed for:

- Women / men
- Overweight / obese
- Light to heavy tendency to venous oedema
- Also for soft connective tissue

Indications²

- Chronic venous disorders: C0s to C4 according to CEAP
- Initial phase after varicose vein therapy
- Oedema during pregnancy
- Postoperative reperfusion oedema
- Post-traumatic, post-operative, occupational, drug-related oedema
- Oedema due to immobility
- Superficial venous thrombosis
- Condition after thrombosis, postthrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Inflammatory dermatoses of the legs
- Sickness, dizziness during pregnancy

Features

- Combines ideal wearing properties and appropriate medical compression
- Your economic solution for a reliable compression therapy



Knee length stocking (AD) with closed toe



Thigh length stocking (AC) with closed toe



Panty (AT) with open toe

¹ Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.
² Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

Product benefits



Clima Fresh



Material

Medium

Compression class

CCL 1, 2

Standard colours

S = Standard



Duomed					
Styles	open toe	AD	AD pte.	AG	AG pte.
		calf	calf petite	thigh	thigh petite
Class 1	Article No.	V140x02-6	V147x02-6	V160x02-6	V150x02-6
Class 2	Article No.	V240x02-6	V247x02-6	V260x02-6	V250x02-6
Styles	open toe	AG /SN	AG pte. / SN	AT	AT pte.
		thigh with topband	thigh petite with topband	panty	panty petite
Class 1	Article No.	V161x02-6	V151x02-6	V110x02-6	V117x02-6
Class 2	Article No.	V261x02-6	V251x02-6	V210x02-6	V217x02-6
		AG	AG pte.		
		thigh with waist attachment right	thigh petite with waist attachment left	thigh with waist attachment right	thigh petite with waist attachment left
Class 1	Article No.	V163x02-6	V162x02-6	V153x02-6	V152x02-6
Class 2	Article No.	V263x02-6	V262x02-6	V253x02-6	V252x02-6
Styles	closed toe	AD	AD pte.	AG	AG pte.
		calf	calf petite	thigh	thigh petite
Class 1	Article No.	V140x12-6	V147x12-6	-	-
Class 2	Article No.	V240x12-6	V247x12-6	-	-
Styles	closed toe	AG /SN	AG pte. / SN	AT	AT pte.
		thigh with topband	thigh petite with topband	panty	panty petite
Class 1	Article No.	V161x12-6		V110x12-6	V117x12-6
Class 2	Article No.	V261x12-6		V210x12-6	V217x12-6

See the last page for the article numbers






Standard styles

Standard sizes	S – XXL
Made -to-measure	not available
Special models	not available
Border	–
Gusset	–
Foot part	With open or closed toe

Material composition


	Elastane (EL)	Polyamide (PA)
CCL 1	37 %	63 %
CCL 2	43 %	57 %

Washing instructions

Wash-and-wear     
 Suitable for washing machine and drier

Measuring Duomed

Size	cB Ankle measurement	cC Calf measurement	cG Thigh measurement	Measurements for silicone topband
S	19 – 21	28 – 34	42 – 57	47 – 54
M	22 – 24	32 – 38	48 – 64	52 – 60
L	25 – 27	36 – 42	54 – 71	56 – 65
XL	28 – 30	40 – 46	60 – 78	62 – 72
XXL	30 – 32	42 – 50	65 – 85	65 – 76



Length

ℓG regular	72 – 83
ℓG petite	62 – 71
ℓD regular	39 – 44
ℓD petite	34 – 38

Toe

For all closed toe stockings:
 only one highly elastic foot size