

Greencare[®]

MEASURED MOBILITY




JUNIOR CHAIR

BRITISH BUILT WHEELCHAIRS

A folding wheelchair for the younger person with in-built stability and strength. The chair can be easily adapted to accommodate growth or changing needs or be easily refurbished for a hospital/care centre environment. The chair is the entry point to access the full range of configurability unique to Greencare and the extensive range of Greencare accessories.

STANDARD FEATURES

- Standard seat widths: 11", 12", 13", 14" on single braced frames for up to 75KG users (double braced frames for added strength can be requested as an option for users to 125KG).
- Standard seat depths: 11", 13", 15", 17" Most frames offer an option to extend seat depth by 2" (50mm) to accommodate growth.
- Seat to ground 17" (option to specify heights 14" to 21").
- Outrigged frame option. (Pic 1) A choice of frames give a straight frame option or an outrigged frame that increases front stability and moves the footrests away from the seat line increasing the space at the front of the chair.
- Rear wheel 24" with other options available, all wheels have a colour coordinated spoke guard for additional safety.
- All junior chairs have a choice of colours for castors, rear tyres and supersoft push rims. Choose from fun pink, myrtle purple, eco green, sea blue, sunflower yellow or pillar box red (Pic 2).
- Height adjustable armrests (Pic 3) allowing adjustment for growth.
- Height adjustable back posts (Pic 4) Not all carers are the same height and therefore the push handles can be height adjusted to meet individual carer's needs.
- A lap belt is fitted as standard. All other posture belts can be fitted on request as can seat positioning aids and power assist options.
- Greencare wheelchairs are crash tested to ISO 7176-19.



IT'S THE LITTLE THINGS THAT MAKE GREENCARE SPECIAL

Chairs are built by Greencare in the United Kingdom and can be customised to meet individual needs. Using quality components designed by the finest British engineers Greencare chairs can be matched by no other as they are the most versatile and uniquely configurable chairs around.

Important dimensions such as seat depths, seat to ground height other than those shown as standard are possible.

You won't be alone when ordering a wheelchair from Greencare as a trained Customer Service Team are on standby to help answer any questions you may have.



Greencare wheelchairs are environmentally friendly, they use the finest aircraft grade aluminium, have a natural anodised finish that minimises unsightly scuff marks commonly experienced on painted surfaces.

We take pride in making long lasting comfortable wheelchairs that are good for the occupant but also work in the environment in which they have to operate, taking into account tight spaces, ramps, small doors and difficult surfaces.

Greencare - a world class product.



ACCESSORIES

Standard ●
Optional ○


SEAT



Lapstrap
With airline type buckle.



Thoracic support
Adjusts for height, width, length on non folding backs.



Extended handles
To aid fitting hub brake levers or other controls.



Extended seat
Brings occupant forward, requires outrigger footrests.



Personalise
With coloured tyres, signwriting or green side panels.



Outrigger frame
Outrigger frame that increases front stability.


WHEELS



Spoke guards
In clear or black plastic.



Greenrim ergonomic design for easy push
Warm to touch and scuff free - 24" only.



Brake lever extension
Fold down lever to aid side transfer.



One arm drive
For 24" wheels please state drive side.



Super heavy duty wheel
24" and 26" wheels tested to 320KG available with hub brakes.



Hub brakes
Available on all wheels sizes. Needs extended handles option.

FRAME




Reinforced back post
Best used when backpost in forward position.



Crutch holder
Left or right fixing.



Oxygen carrier
Uses rear backpost receiver for fixing.



IV pole
1.3m above seat rails. State left or right mount.



Drop frame
Drop frame for ultra-low seat to ground.



Reinforced side panels
Aluminium rigid construction.

ARMREST



Complete tray
With soft edges and cushioned front.




Therapy tray
Non 'handed' cushioned tray with positioning options.



Universal holder
For keeping close handy items, water, TV remote etc.



Tray mount bracket
Unique self centering brackets taking the fuss out of tray mounting.



Extended armrest pad
Extra long arm pad.



Height adjustable armrests

LEGS



ELR with calf pad
Length 23" 580mm. Can be supplied without footplate. Left or right.



Stump supports
Residual limb support with cushioned pad.



Calf straps
Sized to chair.



Adductor pads
Adjustable pad on footrest hanger.



High level release
For users who can't lean forward to release footplate catches.



Dynamic footboard
Carcassone style swinging footplate.


FEET



Angle adjustment footplates



Universal footplates
adjustment 80° to 115° + footplate positioning. Replaces old CAPS 5.



Padded footplates
Soft cushion pad replaces foot plate.


CUSHIONS



Viscoe elastic and high resilience foam core supplied sized to chair in 2" and 3" thickness. The fabric is washable with non slip on base.



Making the right choice of wheelchair is important, that is why we are a preferred supplier to the UK's National Health Service. We have a team on hand to help with the little decisions that will turn a great wheelchair into a perfect choice so if you are choosing a chair for you or someone you care for or love then let us help.



To order call +44 (0)1642 22 33 22
or email info@greencaremobility.com

CARING AND SUPPORTING

Greencare lightweight wheelchairs are amongst the most environmentally friendly wheelchairs in the world and certainly one of the most comfortable, designed by Europe's leading wheelchair designer David Batty and manufactured in the North East England by a dedicated proud workforce.



A VERY SPECIAL WARRANTY

Greencare chairs have a 36 month replacement part warranty on any part found to have a manufacturing defect. Technical manuals, Technical Guide's, User Guide's, assembly drawings, service information, part numbers and more are available on line from the Greencare Website www.greencaremobility.com Servicing and repairs should be conducted by qualified approved repairers.

Greencare in-country support is available in Abu Dhabi, China, Dubai, Eire, Qatar and Saudi Arabia.
© Copyright, FFT Ltd 2015 All rights reserved. Greencare Mobility, Simcox Ct, Middlesbrough TS2 1UU.
Whilst every care was taken in preparing this leaflet E&OE.