

Levamed stabili-tri

Triple-action modular orthosis for the functional care of ankle injuries



Indications

- ankle sprains (grades II to III)
- chronic instability of the talocrural and subtalar joints
- post-operative care after surgery on the ligamentous apparatus of the ankle joint


Mode of action

- modular design supports functional treatment
- adaptable to the three healing phases – for safe and reliable treatment and mobilisation of the patient
- talus lock tongue reduces anterior talar shift
- medial stabilising element gives support and relieves stresses on the lateral side of the foot
- stabilising strap guides the foot in gentle pronation and supports the healing process
- compressive knit improves proprioception and promotes more rapid absorption of haematomas and oedemas

Product benefits





Especially wide opening tongue for comfortable donning and doffing.



Talus lock tongue with cross-lacing offers flexible pressure customisation and simple handling thanks to click function – laces 'click' quickly and reliably into the eyes.



The medial position of the stabilising element protects the injured lateral side of the foot and ensures high comfort in wear.



Clima Comfort



Clima Fresh



DESIGN
AWARD
2015



circumference ankle (cm)	17 – 21	21 – 24	24 – 27	27 – 30
size	I	II	III	IV
article number	K.012.3◊1	K.012.3◊2	K.012.3◊3	K.012.3◊4

◊ = fill in for orientation: 2 = left, 3 = right