M.step®

Functional ankle brace with foam-gel-padding



Indications

- conservative and post-operative care of fibular ligament ankle injuries (grades II to III)
- · chronic instabilities of the ankle
- after distorsions

Mode of action

- rigid shells provide very good stabilisation in pronation and supination movements
- the foam-gel of the padding adapts perfectly to the varying anatomic conditions and provides a uniform compression

Product benefits



Only two straps allow a quick and secure fastening; straps can be taken off the shells and individually fixed.



black





Three length versions available: standard (long version), athletic (short version), kidz (childrens version).



Very light foamgelpadding with soft, skin friendly surface.

| variant | Standard | Athletic | Kidz |
|----------------|-----------|-----------|-----------|
| article number | G.090.009 | G.091.009 | G.093.009 |