

# M.step® Kidz

## Ankle orthosis for stabilisation



### Indications

- conservative and post-operative care of ligamentous injuries of the ankle joint (grades II and III)
- after sprains
- for chronic instability of the ankle joint

### Mode of action

- rigid half shells stabilise pronation and supination movements
- autoadaptive cushion with integrated foam inlays for pressure relief



### Note

From an arm length of more than 29 cm the medi Epico ROMs can be prescribed (see page 103).

size	universal
article number	G.093.009

# medi Walker boot Kidz

## Fixed walker boot



### Indications

- post-operative treatment (rupture of the Achilles tendon, injuries of ligaments, soft-tissues and tendons, fracture of the fibula)
- conservative treatment (fore- and middle foot and tarsus, Achilles tendon, fracture of the fibula)

### Mode of action

- fixes the ankle at 90°
- rocker sole promotes the natural gait
- simple wound inspection and hygiene

Shoe size (EU)	18 – 22	23 – 25	26 – 28	29 – 32
size	S	M	L	XL
article number	G.190.0K1	G.190.0K2	G.190.0K3	G.190.0K4