



## mediven® comfort

mediven comfort is made of wonderfully soft material. The stocking in compression classes 1 and 2 makes the therapy and treatment of venous disorders as comfortable as possible. Thanks to its lanolin finish, it cares for the skin and treats it particularly gently. How do you win over your customers for mediven comfort? Let them feel for themselves!

mediven comfort. The cuddly side of compression.

# mediven<sup>®</sup> comfort

All you feel is softness.  
The cuddly side of compression.

## medi tip<sup>1</sup>:

Due to the material thickness, the fit and the product characteristics, we recommend mediven comfort for:

- Women
- Normal weight / overweight<sup>2</sup>
- Light to severe tendency to venous oedema<sup>2</sup>
- Firm connective tissue

## Indications<sup>3</sup>

- Chronic venous disease stages I to II according to Widmer or C0s – C4 according to CEAP
- Initial phase following varicose therapy
- Congestion problems and oedema during pregnancy
- Postoperative reperfusion oedema
- Post-traumatic, postoperative, occupational, medicinally induced oedema
- Oedema due to immobility
- Superficial venous thrombosis
- Condition after thrombosis, post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Inflammatory dermatoses of the legs
- Nausea and dizziness in pregnancy

## Features

- Particularly soft material for easy donning and doffing
- Added lanolin<sup>4</sup>
- Double-knitted comfort sole for that enhanced feeling of comfort when walking
- Shaping comfort panty top



- Elastic and hard-wearing comfort gusset

## Washing instructions

Easy-care 

Suitable for washing machine and drier.

<sup>1</sup> Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

<sup>2</sup> Use mediven comfort for serious oedema only for patients with normal weight. For people who are overweight and have severe oedema, medi recommends stronger materials to prevent lace-up grooves.

<sup>3</sup> Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

<sup>4</sup> Natural substance used to care dry and chapped skin, has a broad range of applications, for example in child care and baby care.



Comfort panty top



Double-knitted comfort sole



Silicone dot topband

## Product benefits



Clima Comfort



Clima Fresh



Variety



Perfect Fit



## mediven comfort

### open toe

Styles	AD	AD pte.	AG / NoB	AG pte. / NoB	AT	AT pte.
	calf	calf petite	thigh with topband	thigh petite with topband	panty	panty petite
<b>Class 1</b>	Z540x01-7	Z547x01-7	Z569x01-7	Z559x01-7	Z510x01-7	Z517x01-7
<b>Class 2</b>	Z640x01-7	Z647x01-7	Z669x01-7	Z659x01-7	Z610x01-7	Z617x01-7

### closed toe

Styles	AD	AD pte.	AG / NoB	AG pte. / NoB	AT	AT pte.
	calf	calf petite	thigh with topband	thigh petite with topband	panty	panty petite
<b>Class 1</b>	Z540x11-7	Z540x11-7	Z569x11-7	Z559x11-7	Z510x11-7	Z517x11-7
<b>Class 2</b>	Z640x11-7	Z640x11-7	Z669x11-7	Z659x11-7	Z610x11-7	Z617x11-7

See the last page for the article numbers

## Standard styles

Standard	Size I – VII, Length standard and petite, topband standard and wide
Made-to-measure	available
Special models	available, see pages 56 – 73
Border	AD tricot border, narrow silicone dot topband AG Silicone dot topband, Silicone dot topband motif, Plain silicone topband rose, Silicone micro dot topband sensitive All topbands are standard styles.
Gusset	AT Comfort gusset
Foot part	open toe, closed toe, Shoe size 36 – 43

## Material composition

	Elastane (EL)	Polyamide (PA)
<b>CCL 1</b>	32%	68%
<b>CCL 2</b>	35%	65%

## medi tip:



medi Butler donning made easier<sup>1</sup>

medi skin and textile care offer support

More information on the medi accessories from page 112

<sup>1</sup> In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.

## Material

# Medium

## Compression class

# CCL 1, 2

## Standard colours

S = Standard  
M = Made-to-measure

Black

S + M

Anthracite

S + M

Navy

S + M

Cashmere

S + M

Caramel

S + M

Rosé\*

S + M

Bronze\*

S + M

Beige\*

S + M

\* Delivery period up to 10 working days.



Available with sparkling crystals