

mediven[®] for men

It's generally more difficult to convince men to wear medical compression than women. However, the arguments of mediven for men speak for themselves and will help you.

This compression socks looks like a fashionable men's sock and is the perfect accompaniment for men. The stylish sock in a "business ribbed design" feels like a second skin.

mediven for men. For the fashionable business look.

mediven[®] for men

The elegant compression sock.

medi tip¹:

Due to the material thickness, the fit and the product characteristics, we recommend mediven for men for:

- Men
- Normal weight
- No /light tendency to venous oedema
- Also for soft connective tissue





Indications²

- Chronic venous disease stages I and II according to Widmer or COs – C3 according to CEAP
- Postoperative reperfusion oedema
- Post-traumatic, postoperative, occupational, medicinally induced oedema
- Oedema due to immobility
- Superficial venous thrombosis
- Thrombosis prophylaxis in mobile patients
- Inflammatory dermatoses of the legs


Features

- Anatomical toe design with practical right/left markings for a perfect fit and user comfort
- Stylish Business Rib Design
- Added lanolin³
- Clima Fresh for long-lasting freshness on the leg and foot
- Extra-wide and especially soft comfort cuff
- Suitable for shoe sizes 39 – 46

Washing instructions

Easy-care    

Suitable for washing machine and drier.



Knee length sock (AD) in fashionable Business Rib Design with comfort cuffs



Anatomically shaped foot with right and left markings

¹ Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

² Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

³ Natural substance used to care dry and chapped skin, has a broad range of applications, for example in child care and baby care.

Product benefits



Clima Comfort



Clima Fresh



Perfect Fit



cf. Ben Ammar, F. (Dermiscan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

mediven for men

Styles

AD

AD pte.



calf

calf petite

Class	Article No.	3A40x11-7	3A47x11-7
Class 1	Article No.	3A40x11-7	3A47x11-7
Class 2	Article No.	3B40x11-7	3B47x11-7

See the last page for the article numbers

Standard styles

Standard	Size I – VII, Length standard and petite
Special models	not available
Border	AD Tricot border
Foot part	closed toe shoe size 39 – 46

Material composition

	Elastane (EL)	Polyamide (PA)
CCL 1	23%	77%
CCL 2	29%	71%

Material

Light

Compression class

CCL 1, 2

Standard colours

S = Standard

Black
S

Anthracite
S

Navy
S

Brown
S



medi tip:



medi Butler donning made easier¹ medi skin and textile care offer support

More information on the medi accessories from page 112

¹ In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.