



mediven[®] forte

This stocking keeps chronic indications in check. After all, it is all about ensuring that customers with significant vein issues actively participate in life while accepting as few compromises as possible.

The mediven forte delivers maximum pressure stability throughout the day. Its strong material achieves the required depth effect. The stocking is also very hard-wearing and durable.

mediven forte. Pleasant compression has never been so powerful.

mediven® forte

All you feel is protection.
Pleasant compression has never been
so powerful.

medi tip¹:

Due to the material thickness, the fit and the product characteristics, we recommend mediven forte for:

- Women / men
- Overweight / obese
- Light to heavy tendency to venous oedema
- Also for soft connective tissue





Indications²

- Chronic venous disease stages I to III according to Widmer or COs – C5 according to CEAP
- Initial phase following varicose therapy
- Functional venous insufficiency due to obesity
- Congestion problems and oedema during pregnancy
- Postoperative reperfusion oedema
- Post-traumatic, postoperative, occupational, medicinally induced oedema
- Oedema due to immobility
- Superficial venous thrombosis
- Deep leg vein thrombosis
- Condition after thrombosis, post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Inflammatory dermatoses of the legs
- Nausea and dizziness in pregnancy

Features

- Extremely hard-wearing seamless material for optimum wearing comfort
- Maximum pressure stability (high working pressure for a deep functional effect)
- Stable soft cuff at the AD finish
- Opaque
- Clima Fresh over the whole leg for a feeling of freshness around the clock
- Two foot-lengths for series with foot tip:
 - Small foot up to 23 cm
 - Large foot up to 23 cm

Washing instructions

Easy-care    

Suitable for washing machine and drier.



Open toe – detail



Double-knitted fabric



Stable soft cuff at the AD finish

¹ Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

² Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

Product benefits



Clima Fresh

Variety

Perfect Fit



Tested for harmful substances.
www.oeko-tex.com/standard100



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

mediven forte

Styles	AD	AD pte.	AG / NoB	AG pte. / NoB	
	calf	calf petite	thigh	thigh petite	
Class 2	F240x01-7	F247x01-7	F260x01-7	F250x01-7	
Class 3	F340x01-7	F347x01-7	F360x01-7	F350x01-7	
Styles	AG	AG pte.			
	thigh with waist attachment right	thigh with waist attachment left	thigh petite with waist attachment right	thigh petite with waist attachment left	
Class 2	F263x01-7	F262x01-7	F253x01-7	F252x01-7	
Class 3	F363x01-7	F362x01-7	F353x01-7	F352x01-7	
Styles	AT	ATE	AT pte.	AT / U	BT
	panty	one legged panty	panty petite	maternity panty	leggings
Class 2	F210x01-7	made-to-measure	F217x01-7	made-to-measure	made-to-measure
Class 3	F310x01-7	made-to-measure	F317x01-7	made-to-measure	made-to-measure
Styles	CT	ET	AT / H	AT / H pte.	
	capri shorts	cycling short	men's leotard	men's leotard petite	
Class 2	made-to-measure	made-to-measure	made-to-measure	made-to-measure	
Class 3	made-to-measure	made-to-measure	made-to-measure	made-to-measure	

See the last page for the article numbers

Material

Strong

Compression class

CCL 2, 3

Standard colours

S = Standard

M = Made-to-measure

Black

S + M

Anthracite

S + M

Navy

S + M

Cashmere

S + M

Caramel

S + M

Rosé*

S + M

Bronze*

S + M

Beige*

S + M

* Delivery period up to 10 working days.



Available with sparkling crystals

Standard styles

Standard	Size I – VII, Length regular and petite, topband regular and wide	
Made-to-measure	available	
Special model	available, see pages 56 – 73	
Border	AD	Tricot border, narrow silicone dot topband
	AF and AG	Tricot border, Silicone dot topband, Silicone dot topband motif, PI silicone topband rose, Silicone micro dot topband sensitive
	All topbands are standard styles	
Gusset	AT	Sizes I – III groin gusset Sizes IV – VII narrow gusset at the back up to the waist,
	AT / U	with closed panty top
	AT / H	opening horizontal
Foot part	open toe, closed toe	

Material composition

	Elastane (EL)	Polyamide (PA)
CCL 2	21%	79%
CCL 3	26%	74%

medi tip:



medi Butler donning made easier¹

medi skin and textile care offer support

More information on the medi accessories from page 112

¹ In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.