



## mediven plus®

mediven plus has been established for over 40 years and caters for a variety of indications for men and women. The proven medical compression stocking in the compression classes 1, 2 and 3 with additional plus points can be relied upon, no matter what. True to the motto: Reliability that endures.

mediven plus. Compression has never proven so trustworthy.

# mediven plus®

All you experience is reliability.  
Compression has never proven so  
trustworthy.

## medi tip<sup>1</sup>:

Due to the material thickness, the fit and the product characteristics, we recommend mediven plus for:

- Women / men
- Normal weight / overweight
- Light to heavy tendency to venous oedema
- Also for soft connective tissue



## Indications<sup>2</sup>

- Chronic venous disease stages I to III according to Widmer or C0s – C5 according to CEAP
- Initial phase following varicose therapy
- Congestion problems and oedema during pregnancy
- Postoperative reperfusion oedema
- Post-traumatic, postoperative, occupational, medicinally induced oedema
- Oedema due to immobility
- Superficial venous thrombosis
- Deep leg vein thrombosis
- Condition after thrombosis, post-thrombotic syndrome
- Thrombosis prophylaxis in mobile patients
- Inflammatory dermatoses of the legs
- Nausea and dizziness in pregnancy

## Features

- Very hard-wearing and opaque
- Optimal fit thanks to the plus heel, optimised toe, ribbed sole reinforcement and 5 cm wide cuff
- Especially recommended for heavier thighs
- Maxi panty top as standard version

## Washing instructions

Easy-care    

Suitable for washing machine and drier.



Maxi panty top as standard



Plus-heel and sole reinforcement




5-cm-wide cuff

<sup>1</sup> Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

<sup>2</sup> Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.



## Product benefits



Clima Comfort



Clima Fresh



Variety



Perfect Fit



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device – application testing realized under dermatological control, Tunis 2012/2013.

## mediven plus

Styles	AD	AD pte.	AG / NoB	AG / NoB pte.	AT	AT + maxi panty top	AT pte.
	calf	calf petite	thigh with topband	thigh petite with topband	panty	panty with maxi panty top	panty petite
<b>Class 1</b>	4A40x01-7	4A47x01-7	4A61x01-7	4A51x01-7	4A10x01-7	4AZ0x01-7	4A17x01-7
<b>Class 2</b>	4B40x01-7	4B47x01-7	4B61x01-7	4B51x01-7	4B10x01-7	4BZ0x01-7	4B17x01-7
<b>Class 3</b>	4C40x01-7	4C47x01-7	4C61x01-7	4C51x01-7	4C10x01-7	4CZ0x01-7	4C17x01-7
Styles	ATE	ATE pte.	AT / U	AT / U pte.	CT	BT	
	one legged panty	one legged panty petite	maternity panty	maternity panty petite	capri shorts	leggings	
	right	left	right	left			
<b>Class 1</b>	4AE3x01-7	4AE2x01-7	4AV3x01-7	4AV2x01-7	4A20x01-7	4A27x01-7	
<b>Class 2</b>	4BE3x01-7	4BE2x01-7	4BV3x01-7	4BV2x01-7	4B20x01-7	4B27x01-7	
<b>Class 3</b>	4CE3x01-7	4CE2x01-7	4CV3x01-7	4CV2x01-7	4C20x01-7	4C27x01-7	
					made-to-measure	made-to-measure	
Styles	ET	AT / H	AT / H pte.	AG	AG pte.		
	cycling short	men's leotard	men's leotard petite	thigh with waist attachment	thigh with waist attachment petite		
				right	left	right	left
<b>Class 1</b>		4A30x01-7	4A37x01-7	4A63x01-7	4A62x01-7	4A53x01-7	4A52x01-7
<b>Class 2</b>	made-to-measure	4B30x01-7	4B37x01-7	4B63x01-7	4B62x01-7	4B53x01-7	4B52x01-7
<b>Class 3</b>		4C30x01-7	4C37x01-7	4C63x01-7	4C62x01-7	4C53x01-7	4C52x01-7

See the last page for the article numbers

## Material

## Medium

## Compression class

## CCL 1, 2, 3

## Standard colours

S = Standard

M = Made-to-measure

Black

S + M

Anthracite

S + M

Navy

S + M

Cashmere

S + M

Caramel

S + M

Rosé

S + M

Bronze

S + M

Beige

S + M

medi Magenta\*

S + M

White\*

S + M

\* Delivery period up to 10 working days.



Available with sparkling crystals


### Standard styles

Standard	Size I – VII, Length regular and short, topband regular and wide	
Made-to-measure	possible	
Special models	possible, see pages 56 – 73	
Border	AD	tricot border, narrow silicone dot topband
	AF and AG	Silicone dot topband, Silicone dot topband motif, Plain silicone topband rose, Silicone micro dot topband sensitive
	All topbands are standard styles	
Gusset	AT	Size I – III: Crotch gusset Size IV – VII: Narrow gusset, back to waist
	Maxi panty top:	Size I – III: narrow gusset, back to waist Size IV – V: wide gusset, back to waist Size VI – VII: extra-wide gusset, back to waist
	AT / U	with closed panty top
	AT / H	opening horizontal
Foot part	open toe, closed toe, shoe size 36 – 43	

### Material composition

	Elastane (EL)	Polyamide (PA)
<b>CCL 1</b>	25%	75%
<b>CCL 2</b>	31%	69%
<b>CCL 3</b>	37%	63%

### medi tip:



medi Butler donning made easier<sup>1</sup>

medi skin and textile care offer support

More information on the medi accessories from page 112

<sup>1</sup> In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.