Manumed active

Soft wrist support with stabilising splint



Indications

- wrist injuries, e.g. sprains, distortions
- osteoarthritis of the wrist
- instability
- tenosynovitis
- irritative conditions (post-traumatic, post-operative and rheumatic)

Mode of action

- the volar splint stabilises and reduces stress on the wrist
- the circular strap compliments the stress relief provided by the compressive knit and the volar hand splint
- pinpoint pressure relief on the ulnar styloid process enhanced by the anatomically designed strap

Product benefits



The soft knitted end in the lower arm area ensures a pressure-free edge without curling or constricting.



Strap opening for correct positioning of the strap and to relieve strain on the outside of the wrist.



Optimally-sized hand splint offers both large and small hands effective support and stabilisation.



silver











in sizes I + II also for children suitable.

Accessories



The pull-on aid included aids application.



circumference wrist (cm)	13-15	15-17	17-19	19-21
size	1	II	III	IV
article number	K.41 <mark>◊.∆</mark> 01	K.41 <mark>◊.∆</mark> 02	K.41 <mark>◊.</mark> Δ03	K.41 <mark>◊.∆</mark> 04

 \diamond = fill in for orientation: 2 = left, 3 = right/ \triangle = fill in for colour: 2 = silver, 8 = sand