

# Manumed T

Wrist brace with thumb support for immobilisation




## Indications

- distortion of the wrist and / or thumb
- post-operative fixation of the wrist and thumb
- lesions of the collateral ligaments of the basal thumb joint
- tendinitis and tendinopathy of the wrist and / or thumb
- inflammatory phases of rheumatic diseases
- carpometacarpal joint arthrosis

## Mode of action

- fixation and stabilisation of the wrist, the thumb saddle and the basal joint of the thumb


## Product benefits



Easy to fit 'slip-on' system.



Anatomically formed, individually mouldable aluminium splint.



|                                 |           |           |           |           |
|---------------------------------|-----------|-----------|-----------|-----------|
| <b>circumference wrist (cm)</b> | up to 14  | 14 – 18   | 19 – 23   | 24 – 28   |
| <b>size</b>                     | 0         | I         | II        | III       |
| <b>article number</b>           | R.110.◇△0 | R.110.◇△1 | R.110.◇△2 | R.110.◇△3 |

◇ = fill in for orientation: 2 = left, 3 = right / △ = fill in for colour: 6 = grey, 2 = black