

Rhizomed® soft

Splint for immobilisation of
the thumb




Indications

- carpometacarpal arthritis
- ligament injury (e.g. skier's thumb)
- inflammation of the thumb joint

Mode of action

- effective stabilisation of the metacarpal phalangeal and carpometacarpal joint of the thumb
- integrated stabilisation provides reliable immobilisation
- two adjustable thumb straps to fit all hand shapes, providing increased stability
- splint's semi-open shape adapts to various degrees of swelling of the thumb
- thumb-index pinch grip and prehension grasp are maintained

Product benefits



High wearing comfort

The splint has a semi-open shape and is designed in an elastic, breathable fabric, ensuring high wearing comfort and enabling it to be adapted to various degrees of swelling.



Increased stability

The individually adjustable Velcro straps provide additional stability.



Easy handling

Large opening for ease of application and removal, with as little pain as possible.



Reliable immobilisation

The integrated stabilisation provides reliable immobilisation of the injured area.

black



circumference wrist (cm)	13 – 17	18 – 23
size	I	II
article number	R.141.001	R.141.002