



Samsoft 175

New generation




Light, very easy to handle, and small

CE

Samssoft 175

Electric folding!

INNOVATIVE



Disassembly unnecessary
The cross-bar fits in the mast through the remote control.
The cylinder remains protected by the structure.

**USEFUL FOR THE
DISTRIBUTOR OR
THE INSTITUTION**

With its balanced and well-designed structure, it can be **vertically stored**, totally assembled and safely.

A new advantage for a generally bulky product!



Light and easy to handle

Thanks to its minimum size and its strong and light structure, **Samsoft 175** is a manoeuvrable and secure patient lift.

Double-roller wheels,
75-mm Ø.
High-level rolling and
soundness.

Optimisation of the space at
the base foot thanks to a
mechanical system without
play.



Thanks to the ergonomic push-handle, people of all sizes can find their optimum hand position.



POWER SUPPLY

- 1 - Emergency stop
- 2 - Electric emergency lowering
- 3 - Operating indicator



Control box mounted
on the mast



Battery
(2 x 12V / 5,2 Ah / 1,8 kg)



Remote control

All the Samssoft+ straps can be used on the new Samssoft 175.



Comfort polyester strap
Medium size – Ref. SA 639 0200
Large size – Ref. SA 639 0210



BATH netted strap
Medium size – Ref. SA 639 0100
Large size – Ref. SA 639 0110



NETTED HAMMOCK strap
One size - Ref. SA 639 0500



TOILET NETTED HAMMOCK strap
One size - Ref. SA 639 0610



U-shaped strap with headrest
One size - Ref. SA 639 0800


SPECIFICATIONS

All the photographs are not contractual, and some components may change for manufacturing or improving purposes.

Reference for Samssoft 175:	SR2810000
Lifting capacity:	175 kg
Maximum lifting height:	157 cm
Minimum floor/base height for the hooks:	55 cm
Overall width of the base:	61 cm
Overall length:	128 cm
Base height:	9 cm
Inner feet spacing:	51 cm to 94 cm
Patient lift weight:	36 kg



**On structure,
cylinder,
and control box**

Complies with NF - EN-standard / Cylinders IP X4 - 

DUPONT MEDICAL

Parc d'activités Eiffel Energie - Site de Ban la Dame
BP 19 - 54390 FROUARD (France)
Tel. +33 (0)3 83 495 495 - Fax. +33 (0)3 83 495 496
<http://www.dupont-medical.com>

SAMERY®
A DUPONT MEDICAL trademark