
Spinomed[®]

Back orthosis for vertebral extension in osteoporosis and strengthening of the core musculature



Indications

- osteoporotic bone collapse in the thoracic and lumbar spine
- juvenile Scheuermann's disease
- hyperkyphosis with chronic back pain

Mode of action

- belt system and back support cause activation, by triggering biofeedback
- reduction of kyphosis from osteoporosis
- increase in mobility by enhancing muscle strength
- reduction of pain





Ergonomically shaped shoulder straps make it easy to put on and increase wearing comfort.

grey



The soft material with inset holes increases air circulation and is comfortable to wear.



The new strap loops help prevent twisting of the straps.

Studies



The effect of Spinomed back orthosis has been confirmed by randomised Studies. Copies available upon request (article 99E22).



back length (cm)	up to 43	44 – 49	50 – 55	56-61	62 – 67
size	XS	S	М	L	XL
article number	5.669.051	5.669.052	5.669.053	5.669.054	5.669.055