

Sliding Equipment



promoting safer patient handling



Lateral Transfer Slide

The Lateral Transfer Slide has a smooth surface and non-slip base allowing a patient to be moved from one surface to another without the need to be lifted. Used in conjunction with the Transfer Sheet, correct posture during the transfer is maintained. The Lateral Transfer Slide has tapered edges to allow for easy insertion and removal.

LOCO-119 (185cm x 65cm) Lateral Transfer Slide

Transfer Sheet

A full body-length sheet with six long handholds which can be used with the Maxi Mover or a bridging board to aid sliding transfers between two surfaces. The long handholds eliminate the need to bend and stretch to effect the transfer. The Transfer Sheet can also be used in an emergency situation when a patient has fallen on the floor where a hoist is unavailable or inappropriate. Weight tested to 127kg/20 stones.

LOCO-009 (190cm x 100cm) Transfer Sheet

Deluxe Transfer Sheet

Used in the same way as the Transfer Sheet above, the Deluxe Transfer Sheet incorporates a series of extra handholds which allows for more precise positioning of a patient during a transfer. The extra handholds also increase control during the handling task by enabling staff to maintain correct posture during the transfer.

LOCO-205 (190cm x 100cm) Deluxe Transfer Sheet

Quilted Slide Sheets

Available in a variety of sizes and used for turning and repositioning patients in bed, for transfer between hospital bed/ trolley/examination couch or operating table. Supplied complete with a waterproof cover, which can be used independently as a slide sheet.

LOCO-004	(190cm x 65cm)	Maxi Mover
LOCO-013	(195cm x 70cm)	Maxi Cover
LOCO-066	(190cm x 130cm)	Maxi Mover Plus
LOCO-067	(195cm x 140cm)	Maxi Plus Cover
LOCO-005	(70cm x 65cm)	Mini Mover
LOCO-012	(72cm x 70cm)	Compact Cover



