

medi Epico ROMs

Elbow brace for mobilisation with immobilisation of the PRU joint




Indications

- treatment of elbow luxations (conservative / post-operative)
- orthotic care for medial or lateral epicondyle fractures
- treatment of fractures of the radial head (conservative / post-operative)
- treatment of distal humeral fractures (conservative / post-operative)
- post-operative care after prosthesis implantation

Mode of action

- orthosis for individual early functional treatment
- movement is limited to 0° – 120° in the hinge, to guarantee movement flexibility appropriate to the stage by early mobilisation – this enhances joint mobility and accelerates healing
- flexion / extension restricted to the range 0° – 120° – adjustable in steps of 10°

Product benefits



Limitation 0° - 120° serves to immobilise the elbow joint by fixation of areas of the upper and lower arm.

Flat Quick Set hinge pleasant to wear, even under clothing. The hinge can be simply adjusted by pressing a button.

Padding made of innovative foam and fits tightly and reliably – the fluted surface slightly compresses the muscle receptors, stimulating metabolism and preventing muscle tension.

Aluminium construction can be precisely shaped and adapted to the shape of any arm.

Length of the hand cushion can be adjusted – the orthosis is securely attached to the arm.

grey / graphit



medi Kidz – also available as Kidz-Version.



length of forearm (cm)	up to 29	29 – 42
size	S	L
article number	R.032.◇02	R.032.◇04

◇ = fill in for orientation: 2 = left, 3 = right