

medi Humeral fracture brace

Humeral fracture brace



Indications

- post-operative care of subcapital or shaft fractures to the humerus
- conservative care of humeral fractures

Mode of action


- stabilisation of fracture due to metal bars and strap system
- slight compression
- easy to apply and remove for early rehabilitation

Product benefits

 black



Rigid bars removable for custom fit and x-ray.



Fully adjustable upper arm belt


Option to reposition or remove upper arm belts for flexible range of application and compression.



Fully adjustable upper arm belt

Option to reposition or remove upper arm belts for flexible range of application and compression.

Extension pad for larger circumference measures available as spare part
article number
R.E90.039.



circumference upper arm (cm)	21–28	26–34	31–40
size	S	M	L
article number	R.091.◇12	R.091.◇13	R.091.◇14

◇ = fill in for orientation: 2 = left, 3 = right