

protect.Ankle foot orthosis

Ankle orthosis with adjustable strap system

Indications

- foot drag due to hemiplegia
- drop foot due to peroneal paralysis

Mode of action

- the Ankle foot orthosis raises your foot while the leg swings forward and prevents your foot from tipping down



Product benefits

anatomically preshaped

lightweight

waterproof

easy to clean

can be worn with a shoe

ice may be applied with the orthosis in place

protect.Ankle foot orthosis (size of shoe)

| size | shoe size EU | article number left | article number right |
|------|-----------------|------------------------|-------------------------|
| S | 35 – 37 | F.AFO.202 | F.AFO.302 |
| M | 37 – 39 | F.AFO.203 | F.AFO.303 |
| L | 39 – 41 | F.AFO.204 | F.AFO.304 |
| XL | 41 – 44 | F.AFO.205 | F.AFO.305 |