

# protect.CAT Walker

Comfortable walker for ankle and foot injuries – in a lightweight shell construction.

## Indications

- postoperative: achilles' tendon rupture, post ligament reconstruction, treatment of soft tissue injuries and tendon injuries
- conservative: achilles' tendon injuries (2nd and 3rd grade), ligament lesions, tenosynovitis
- stable foot and joint fractures: metatarsal fracture, talus fracture, malleolar fracture
- after ligament repair

## Mode of action:

- safe and effective immobilisation
- four strap system
- ensures high stability and reduces anterior talar shift
- available as optional extra: heel wedge sets for Achilles' tendon rupture



## Product benefits



air system with medial and lateral air passage around the ankle and up to the forefoot



air pump and ventilation valve for easy handling



easy adjustment and removal of toe protector



flat, rounded sole for natural, stable gait

## Accessories



Hygiene cover prevents the spread of dirt, dust and other contaminants or germs, which can occur when using the walker both outside and indoors. Available for all medi walker products (article number G.E19.009).

**new article numbers**

**protect.CAT Walker** (size of shoe)

size	EU	article number
S	36 – 40	G.950.012
M	40 – 44	G.950.013
L	44 – 47	G.950.014