# protect.CAT Walker short

Comfortable walker for ankle and foot injuries – in a lightweight shell construction.

## Indications:

- postoperative: post ligament reconstruction and/or refixation, for treatment of soft tissue injuries and tendon injuries
- conservative: ligament lesions 2nd and 3rd grade, tenosynovitis
- stable foot and joint fractures: metatarsal fracture, talus fracture, malleolar fracture

## Mode of action:

- safe and effective immobilisation
- three strap system ensures high stability and reduces anterior talar shift
- available as optional extra: heel wedge sets for Achilles' tendon rupture



#### Product benefits



air system with medial and lateral air passage around the ankle and up to the forefoot



air pump and ventilation valve for easy handling



easy adjustment and removal of toe protector



flat, rounded sole for natural, stable gait



# protect.CAT Walker short (size of shoe)

size	EU	article number
S	36 – 40	G.950.112
M	40 – 44	G.950.113
L	44 – 47	G.950.114

### Accessories



Hygiene cover prevents the spread of dirt, dust and other contaminants or germs, which can occur when using the walker both outside and indoors.

Available for all medi walker products (article number G.E19.009).