

protect.Dorsofix

Semi-rigid thoracolumbar corset

Indications

- acute and chronic middle and lower back pain
- vertebral fracture recovery period
- lumbar radiculitis
- herniated discs
- vertebral arthrosis
- curvature of the thoracic spine
- mild spondylolisthesis
- post-operative period after lumbar discectomy surgery
- recovery period following spinal stabilisation operations

Mode of action

- adjusts perfectly to the body's anatomical shapes, providing optimum stability of the thoracolumbar segment



Product benefits



made in an innovative material combining the characteristics of cotton and Lycra




equipped with shoulder straps, bars, Velcro closures, elastic straps and a front polyethylene reinforcement

pleasant to wear due to breathable material

easy adjustment to the thoracic spine due to the aluminium bars

available as men's and ladies' version

protect.Dorsofix



measurement
circumference hip

| size | cm | article number women |
|------|-----------|----------------------|
| XS | 78 – 86 | K.660.821 |
| S | 86 – 96 | K.660.822 |
| M | 96 – 106 | K.660.823 |
| L | 106 – 116 | K.660.824 |
| XL | 116 – 126 | K.660.825 |
| XXL | 126 – 136 | K.660.826 |
| size | cm | article number men |
| XS | 71 – 81 | K.660.811 |
| S | 81 – 91 | K.660.812 |
| M | 91 – 101 | K.660.813 |
| L | 101 – 111 | K.660.814 |
| XL | 111 – 121 | K.660.815 |
| XXL | 121 – 131 | K.660.816 |