

# protect.Knee immobilizer universal

## Knee immobilisation brace at 0°

### Indications

- post-operative immobilisation
- pre-operative pain conditions
- following patella dislocations

### Mode of action

- immobilisation of the knee at 0°



### Product benefits

individual, rapid adjustment to fit varying leg circumferences

removable medial, lateral and posterior stabilisation rods

comfortable, washable, skin friendly material

### protect.Knee immobilizer universal

sizes	length 40cm	length 50cm	length 60cm
universal	5.652.S49	5.652.S59	5.652.S69