

# protect.ROM telescopic

Universal knee splint with flexion / extension limitation and telescopic length adjustment

## Indications

- post-operative immobilisation
- for secure stabilisation of injured or operated collateral and / or cruciate ligaments

## Mode of action

- collateral stabilisation of the knee joint
- post-operative immobilisation
- controlled mobilisation
- limitation in steps of 10° to 15°
- flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
- extension: 0°, 10°, 20°, 30°
- immobilisation: 0°, 10°, 20°, 30°



## Product benefits



telescopic system: bar length adjustable from 42 – 52 cm

patented Quick-Set hinge

fits up to 80 cm thigh circumference

## protect.ROM telescopic

size	length (cm)	article number
universal	42 – 52	G.185.009