

# protect.ROM / protect.ROM cool

Universal knee splint with flexion- / extension limitation and options for immobilisation

## Indications

- post-operative immobilisation
- for secure stabilisation of injured or post-op collateral and / or cruciate ligaments

## Mode of action

- collateral stabilisation of the knee joint
- post-operative immobilisation
- controlled mobilisation
- immobilisation in steps of 10° to 15°
- flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
- extension: 0°, 10°, 20°, 30°
- immobilisation: 0°, 10°, 20°, 30°



protect.ROM



protect.ROM cool

## Product benefits



easy to use patented Quick-Set hinge

option to shorten the side bars due to integrated break lines

skin friendly padding which can be adjusted according to patient needs by cutting

protect.ROM		
version	length (cm)	article number
protect.ROM	63	G.610.509
protect.ROM short	57	G.610.519
protect.ROM cool	63	G.620.509
protect.ROM cool short	57	G.620.519