

protect.Walker boot

Fixed Walker

Indications

- foot or toe injuries or surgery
- grade I – III ankle sprains
- stable or internally fixed fractures of the foot or ankle

Mode of action

- rocker sole encourages natural gait



Product benefits



also available as short version

available in 4 sizes, adjustable circumference

velcro system facilitates wound inspection and personal hygiene

fixed Walker has choice of standard or short length

Accessories



Hygiene cover prevents the spread of dirt, dust and other contaminants or germs, which can occur when using the walker both outside and indoors. Available for all medi walker products (article number G.E19.009).

protect.Walker boot (size of shoe)

size	EU	article number
S	36 – 39	G.920.012
M	40 – 44	G.920.013
L	44 – 46,5	G.920.014
XL	46,5 +	G.920.015

protect.Walker boot short (size of shoe)

size	EU	article number
S	36 – 39	G.920.022
M	40 – 44	G.920.023
L	44 – 46,5	G.920.024
XL	46,5 +	G.920.025

also available:

protect.Air booties and protect.Walker wedges (size of shoe)

size	EU	article number protect.Air booties	article number protect.Walker wedges
S	36 – 39	G.E12.012	G.E11.132
M	40 – 44	G.E12.013	G.E11.133
L	44 – 46,5	G.E12.014	G.E11.134
XL	46,5 +	G.E12.015	G.E10.135