

# IDSoft Evolution

Accurate conception of comfort and positioning for various pathologies

 $\epsilon$ 



# The right identity of positioning

## Two versions

**IDSoft** Evolution



#### **IDSoft** Evolution Prime



# and options



Long wings headrest



standard shelf



Viscoelastic foam seat cushion



Swiwelling cradled armrest



Antitip castors



Trunk support and abduction wedge

## and comfort

Thanks to the frame design, the wheelchair is steady, whatever the backrest and seat angles are (through gas cylinder).





A very functional multiple-adjustment legrest

Large, adjustable in heigh, width and depth calf cushion.

Rotation-, and front-to-back adjustable board.







Various adjustments, maximum adaptability for comfort, and user position and security.



Width between both armrests (= seat width)

Continuous adjustment:

On T 39 frame from 39 to 44 cm
On T 44 frame from 44 to 49 cm
On T 49 frame from 49 to 54 cm



Long and curved armrest, heightadjustable on 70 mm.

Whole detachable armrest through a push-button.



Continuous depth-adjustable seat, through sliding, on 90 mm (42 to 51 cm).



**Seat reclining** from –2° (Trendelenburg) to +20° with large rear wheels, and from -2° to +25° with small rear



Heigh adjustable backrest (on 10 cm).



Reclining backrest from 90° (vertically) to 120°.

#### Floor/seat plate heigh available:

Floor/seat plate heigh	Rear Wheels	Rear wheel position	Fork position	Front wheel
45 cm	24"	2	2	8"
42.5cm	24"	1	2	6"
42.5cm	22"	2	2	6"
40 cm	22"	1	1	6"
45 cm	12,5"	5	2	8"
40 cm	12,5"	4	1	6"

### COMPACT SIZE

## MANOEUVRABILITY

#### **STABILITY**

Overall width	Overall length without legrests	Overall length with legrests	Wheelchair weight	Maxi. user's weight
Frame width + 20 cm T39 at 44 = 59 cm T44 at 49 = 64 cm T49 at 54 = 69 cm	82 cm	108 cm	From 34,2 kg	140 kg

The photographs are not contractual, and some components may change for manufacturing or improving purposes.

 $\epsilon$ 



BP 19 - Parc d'activités Eiffel Energie Site de Ban la Dame 54390 FROUARD (France) Tel. +33 (0)3 83 495 495 - Fax. +33 (0)3 83 495 496 www.dupont-medical.com

