

mediven elegance[®]

Is this really a medical compression stocking? For sure! This is a model that emphasises a woman's beauty. The mediven elegance pampers the legs with its compressive effect. In doing so, it nestles against the leg like a second skin. It can be combined with any item of clothing and perfectly accentuates the female figure. True to the motto: Beauty that can be seen.

mediven elegance. Effective compression has never looked so good.

mediven elegance®

All you see is beauty.
Effective compression has never
looked so good.

medi tip¹:

Due to the material thickness, the fit
and the product characteristics, we
recommend mediven elegance for:

- Women
- Normal weight
- No tendency to venous oedema
- Firm connective tissue






Indications²

- Chronic venous disease stage I
according to Widmer or COs – C2
according to CEAP
- Initial phase following varicose therapy
- Superficial venous thrombosis
- Thrombosis prophylaxis in mobile
patients
- Inflammatory dermatoses of the legs
- Nausea and dizziness in pregnancy

Features

- Transparent fine stocking style
- Soft Elastic foot tip guarantees more
freedom for the toes
- Wide and soft cuff for knee length style
- Smaller heel for wrinkle-free fit
- Attractive topband motifs for
high-length stockings

Washing instructions

Easy care     

Suitable for washing machine and drier.



Soft elastic foot tip



Knee length stocking (AD)
with soft cuff



Plain silicone topband rose (FBR)

¹ Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

² Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

Product benefits



cf. Ben Ammar, F. (DermScan): Clinical evaluation of the skin compatibility of a medical device - application testing realized under dermatological control, Tunis 2012/2013.

mediven elegance

Styles	AD	AD pte.	AG	AG pte.	AG / FBR
	calf narrow silicone topband	calf petite	thigh	thigh petite	thigh with topband
Class 1	0140x11-7	0147x11-7	0160x11-7	0150x11-7	016Hx11-7
Class 2	0240x11-7	0247x11-7	0260x11-7	0250x11-7	026Hx11-7
Styles	AG pte. / FBR	AT	AT pte.	AT / U	AT / U pte.
	thigh petite with topband	panty	panty petite	maternity panty	maternity panty petite
Class 1	015Hx11-7	0110x11-7	0117x11-7	0120x11-7	0127x11-7
Class 2	025Hx11-7	0210x11-7	0217x11-7	0220x11-7	0227x11-7

See the last page for the article numbers

Standard styles

Standard	Size I – VII, Length regular and petite, topband regular and wide
Made-to-measure	available
Special models	available, see pages 56 – 73
Border	AD tricot border, narrow silicone dot topband
	AF and AG tricot border, Silicone dot topband, Silicone dot topband motif, Plain silicone topband rose, Silicone micro dot topband sensitive
	All topbands are standard styles.
Gusset	AT Comfort gusset
	AT / U with closed panty top
Foot part	open toe, closed toe, Shoe size 36 – 43

Material composition

	Elastane (EL)	Polyamide (PA)
CCL 1	34%	66%
CCL 2	39%	61%

medi tip:



medi Butler donning made easier¹

medi skin and textile care offer support

More information on the medi accessories from page 112

¹ In general, the following applies for donning and doffing aids: Rabe E et al. Medical compression therapy of the extremities with medical compression stockings (MCS), phlebological compression bandages (PCB), and medical adaptive compression systems (MAC): S2k guideline of the German Phlebology Society (DGP) in cooperation with the following professional associations: DDG, DGA, DGG, GDL, DGL, BVP. Hautarzt 2021;72(2):137-152.

Material

Light

Compression class

CCL 1, 2

Standard colours

S = Standard
M = Made-to-measure

Black

S + M

Anthracite

S + M

Navy

S + M

Cashmere

S + M

Caramel

S + M

Rosé

S + M

Bronze

S + M

Beige

S + M


medi Magenta*

S + M

White*

S + M

* Delivery period up to 10 working days.



Available with sparkling crystals