



mediven sheer & soft®

medi's most translucent compression stocking. Try a major advance in compression therapy. mediven sheer & soft takes comfort to a whole new level. Yet these compression stockings are so sheer that only you will know that they are medical compression stockings.

mediven sheer & soft®

medi's most translucent compression stocking.

medi tip¹:

Due to the material thickness, the fit and the product characteristics, we recommend mediven sheer & soft for:

- Women
- Normal weight
- No tendency to venous oedema
- Firm connective tissue

Indications²

- Chronic venous disorders: C0s to C2 according to CEAP
- Initial phase after varicose vein therapy
- Superficial venous thrombosis
Thrombosis prophylaxis in mobile patients

Features


- Extremely sheer, not distinguishable from fashion stockings
- Soft and gentle on the skin
- Available with open toe



Attractive topband



Knee length stocking (AD)



Panty (AT)



Open toe – detail

¹ Please note that this recommendation does not replace a doctor's diagnosis or professional advice from a specialist medical retailer.

² Selecting the ideal compression stockings depends on the indication and the individual needs of the patients.

Product benefits



Material


Light

Compression class

15 – 20 mmHg,
20 – 30 mmHg

Standard colours

S = Standard



mediven sheer & soft

Styles	AD	AD pte.	AG	AG pte.
	calf	calf petite	thigh	thigh petite
15 – 20 mmHg	O140x11-7	O147x11-7	O161x11-7	O151x11-7
20 – 30 mmHg	O240x11-7	O247x11-7	O261x11-7	O251x11-7
Styles	AT	AT pte.		
	panty	panty petite		
15 – 20 mmHg	O110x11-7	O117x11-7		
20 – 30 mmHg	O210x11-7	O217x11-7		

See the last page for the article numbers

Standard styles

Shoe size 35 – 42

Material composition

	Polyamide (PA)	Elastane (EL)
15-20 mmHg	64 %	36 %
20-30 mmHg	64 %	36 %

Washing instructions

Wash-and-wear Suitable for washing machine and drier

Measuring mediven sheer & soft

Size	Circumference in cm						
	I	II	III	IV	V	VI	VII
cG	42 – 57	45 – 61,5	49 – 66,5	52 – 73	56 – 79	60 – 84,5	64 – 87,5
cD	26 – 33,5	28 – 37	30 – 40	32 – 43	33,5 – 46,5	35,5 – 49,5	37,5 – 52,5
cC	28 – 35,5	30 – 39	32,5 – 42	35 – 45	37 – 48,5	39,5 – 51,5	41,5 – 54,5
cB	18 – 20	20,5 – 22,5	23 – 25	25,5 – 27,5	28 – 30	30,5 – 33	33,5 – 36
Length							
Calf stocking / panty petite:			34 – 38 cm		Thigh stocking / panty petite:		62 – 71 cm
Calf stocking / panty standard:			39 – 44 cm		Thigh stocking / panty standard:		72 – 83 cm