

Locomotion





Unit Width:

Unit Length:

NEW

The Locomotion facilitates a dignified transfer enabling a person to travel safely over short distances from one seated position to the next.

Promotes safe standing

Rear wheel steer for increased movability in confined spaces

Allows single carer operation*

Single handed shin pad height adjustment - with independent swivel and lateral movement

Single point brake


Central mast for easy leg positioning

Dual height accessory attachment points

200kg SWL (safe working load)

* Subject to individual assessment.

LOCO-501 The Locomotion



Used to assist the user to a standing position and to support once stood.

Small (red)

Medium (yellow)

Large (green)

Patent: GB1908712.1

LOCO-514 Locomotion Assist Belt - Small (Waiste 23" - 48")

LOCO-515 Locomotion Assist Belt - Medium (Waiste 33" - 58")

LOCO-516 Locomotion Assist Belt - Large (Waiste 43" - 78")

Rota Stand Compact

The Rota Stand Compact combines all the benefits of the Rota Stand and the Rota Stand Solo in one product with the added value of being lightweight (only 10kg) and affordable.

Foldable compact design reduces storage and transport issues

Lightweight: easy to manoeuvre and lift in and out of vehicles

Height adjustable frame and shin pad to suit all clients

Central wheels allow easy positioning

Stabilizing bar allows a one person transfer

200kg weight limit

The folding and lightweight nature of the Rota Stand Compact is ideal for community assessments by therapists who often have to lift equipment in and out of their vehicles.

Equally, a more independent lifestyle is encouraged for clients who wish to travel to friends and relatives. The Rota Stand Compact fits easily into the boot of a car.

LOCO-600 Rota Stand Compact








NEW



114cm



Rota Stands



promoting safer patient handling

Rota Stand



Allows a safe, smooth and controlled swivel transfer



Reduces the need to "hold" a patient during a transfer



Encourages an independent and dignified transfer

The Rota Stand has been developed to allow an assisted transfer from one seated position to another. The frame, shin pad and turn table combine to give optimum safety to both carer and patient during the transfer.

- Height adjustable frame to suit all
- Foot pedal break to secure turn-table
- Central wheels for easy positioning
- Angled shin pad for correct standing position
- 160kg weight limit
- Height adjustable shin pad

The central wheels are a unique innovation, because they allow the Rota Stand to be rolled up to the patient or client even after their feet have been positioned on the product.



Patent: GB2389099B

LOCO-102

Rota Stand

Rota Stand Solo

Height adjustable

shin pad



Stabilizing bar to allow a one person transfer



Height adjustable frame to suit all clients



Central wheels for easy positioning



Large handlebar for a secure grip



200kg weight capacity

The Rota Stand Solo revolutionises sit to stand and pivot transfers. The patented T Bar mechanism in the base (see picture) allows a sit to stand and pivot transfer to be carried out by a single nurse or carer thus reducing the time it takes to carry out the task as well as halving the cost of having two carers/nurses present during the transfer.

- The nurse/carer does not need to counterbalance the Rota Stand Solo as the client pulls themselves to a standing position.
- The nurse/carer can now assist the client to stand before rotating the client to a new seated position.



LOCO-311 Rota Stand Solo