Stabimed[®]

Short functional soft brace with flexion / extension limitation



Indications

- · collateral lesions and instabilities
- mild to moderate instabilities of the knee joint
- chronic instability
- meniscus injuries
- · structural relaxation
- traumatic or degenerative deformities of the knee joint

Mode of action

- collateral stabilisation of the knee joint due to side bars
- prevention of hyperextension and limitation of range of motion due to adjustable physioglide hinges
- extension limitation at: 0°, 10°, 20°, 30°, 45°
- flexion limitation at:
 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- immobilisation at:
 0°, 10°, 20°, 30°, 45°
- · compressive and proprioceptive effect





Very comfortable wearing due to persistent net material on the back.



grey





Tool-free wedge replacement.



wrap version



The physioglide hinge mimics the rolling / gliding movement of the knee and ensures a secure fit.



circumference thigh (cm)	31-37	37 – 45	45 – 51	51 – 57	57 – 65	65 – 74	74 – 84
size	XS	S	М	L	XL	XXL	XXXL
article number	G.070.041	G.070.042	G.070.043	G.070.044	G.070.045	G.070.046	G.070.047

Following wedges are pre-installed on the brace at the point of delivery: extension 0° / flexion 90°.