

Indications

• epicondylitis (e.g. tennis or golfer's elbow)

Mode of action

 reduction of the stress on the inflamed muscle and tendon attachments in the elbow due to the direct press from the pad onto the lower arm musculature

Product benefits



Flexible pad: can be used anywhere, on the left and right, as well as for golfer's and tennis elbow.



grey/grey



magenta/grey



lime/grey



New wave structure and pad with 3D profile improve the removal of lymph fluid.



Soft edge: gentle on the skin, does not dig in.



Precise pressure adjustment thanks to the numbering scale on the strap.



Even easier to put on thanks to the support strap.

Usable until 36 cm circumference lower arm

size universal
article number R.110.01∆

 Δ = fill in for colour: 0 = grey / grey, 1 = magenta / grey, 2 = lime / grey