

ET800



TRACION-u Traction Unit



Features Description

• FEATURES

Tracion-u® ET800 is the result of thousands of feedbacks from tens of thousands of physiotherapists over the past decade. The Retro look of Traç Pod® is also one of the many features that physiotherapists requested to retain for good. Instead of a touch screen LCD panel and working thru complicated windows behind windows, physiotherapists recommended the easy-to-adjust "HOLD" and "REST" knobs as they have been proven to be one of the fastest and effective way of adjustment of traction force during treatment. The "HOLD" and "REST" TIME can be adjusted in real-time as well as the FORCE. Once set, press START and the new setting will take effect in the next coming cycle. The unit is equipped with a "Patient Switch" for Full Control of an Emergency Situation, all at the hands of a patient. The Tracion-u® ET800 works best on Traçion® EU200 however it also can be fitted on tables of other make

Controls Behind



For Cervical Traction Treatment, CT850 is available as an optional add-on part. To enable this function, a switch behind the unit can be located to switch between (1) Lumbar Traction and (2) Cervical Traction. Weight reading of kgs and lbs can also be switched via a switch behind the unit.

With the ET800, CT850 and EU200, you have the ultimate traction system in your home, clinic, hospital or rehab center.

Other Optional Accessories



Product Name Order No
HEAD HALTER HH-1



Product Name Order No
THORACIC TRACTION BELT TB-1



Product Name Order No
LUMBAR TRACTION BELT LB-1



Product Name Order No
SPREADER BAR SB-1



Product Name Order No
COUNTER TRACTION STRAPS (PAIR) CS-1



Product Name Order No
CORD EXTENSION CX-1



Product Name Order No
FLEXION STOOL FS-1



• 1 Year Full Manufacturer Warranty